SphygmoCor e malattia renale

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Bibliografia

- 1 Kidney and Urological Disease Statistics for the United States. National Kidney and Urological Diseases Information Clearinghouse, February 2004, NIH Publication No. 04-3895.
- 2 De Nicola L, Minutolo R, Chiodini P, *et al.* Global approach to cardiovascular risk in chronic kidney disease: reality and opportunities for intervention. *Kidney Int* 2006;69:538-545.

- 3 Sarnak MJ, Levey AS, Schoolwerth AC, *et al.* Kidney disease as a risk factor for development of cardiovascular disease. A Statement from the American Heart Association Councils on Kidney in Cardiovascular Disease, High Blood Pressure Research, Clinical Cardiology, and Epidemiology and Prevention. *Circulation* 2003;108:2154-2169.
- 4 U.S. Renal Data System, USRDS 2004 Annual Data Report: Atlas of End-stage renal disease in the United States. National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases, Bethesda, MD, 2004.
- 5 American Heart Association. Heart Disease and Stroke Statistics – 2004 Update. Dallas, Tex.: American Heart Association, 2003.
- 6 Blacher J, Guerin AP, Pannier B, *et al.* Impact of aortic stiffness on survival in end-stage renal disease. *Circulation* 1999;99:2434-2439.
- 7 London GM, Blacher J, Pannier B, *et al.* Arterial wave reflections and survival in end-stage renal failure. *Hypertension* 2001;38:434-38.
- 8 Guerin AP, London GM, Marchais SJ, *et al.* Arterial stiffening and vascular calcifications in end-stage renal disease. *Nephrol Dial Transplant* 2000;15:1014-1021.
- 9 Haydar AA, Covic A, Colhoun H, Rubens M, Goldsmith DJA Coronary artery calcification and aortic pulse wave velocity in chronic kidney disease patients. *Kidney Int* 2004;65:1790-1794.
- 10 Covic A, Gusbeth-Tatomir P, Goldsmith DJA. Arterial stiffness in renal patients: An update. *American Journal of Kidney Diseases* 2005;45:965-977
- 11 Covic A, Mardare N, Gusbeth-Tatomir P, *et al.* Increased arterial stiffness in children on haemodialysis. *Nephrol Dial Transplant* 2005; October 12.
- 12 Liao D, Arnett DK, Tyroler HA, *et al.* Arterial stiffness and the development of hypertension. The ARIC study. *Hypertension* 1999;34:201-206.
- 13 Laurent S, Boutouyrie P, Asmar R, *et al.* Aortic stiffness is an independent predictor of all-cause and cardiovascular mortality in hypertensive patients. *Hypertension* 2001;37:1236-1241.
- 14 Brooks B, Molyneaux L, Yue DK. Augmentation of central arterial pressure in Type 1 diabetes. *Diabetes Care* 1999;22:1722-1727.
- 15 Wilkinson IB, MacCallum H, Rooijmans DF, *et al.* Increased augmentation index and systolic stress in Type 1 diabetes mellitus. *QJM* 2000;93:441-8.
- 16 Shram MT, Henry R, van Dijk R, *et al.* Increased arterial stiffness is impaired in glucose metabolism and Type 2 diabetes. The HOORN study. *Hypertension* 2003;43:176-181.
- 17 Cruickshank K, Riste L, Anderson SG, *et al.* Aortic pulse-wave velocity and it's relationship to mortality in diabetes and glucose intolerance. *Circulation* 2002;106:2085-2090.
- 18 McLeod A, Uren AL, Wilkinson AB, *et al.* Non-invasive measures of pulse wave velocity correlate with coronary arterial plaque load in humans. *J Hypertens* 2004;22:363-368.
- 19 Wilkinson IB, Prasad K, Hall IR, *et al.* Increased central pulse pressure and augmentation index in subjects with hypercholesterolemia. *J Am Coll Cardiol* 2002;39:1005-1011.
- 20 Mourad JJ, Pannier B, Blacher J, *et al.* Creatinine clearance, pulse wave velocity, carotid compliance and essential hypertension. *Kidney Int* 2001;59:1834-41.

- 21 Briet M, Bozec E, Laurent S, *et al.* Arterial stiffness and enlargement in mild-to-moderate chronic kidney disease. *Kidney Int* 2006;69:350-7.
- 22 Lacy P, Carr SJ, O'Brien MB, *et al.* Reduced glomerular filtration rate in pre-dialysis chronic kidney disease in patients is associated with impaired baroreceptor sensitivity and reduced vascular compliance. *Clin Sci* 2006;110:101-108.
- 23 London GM, Pannier B, Vicaud E, *et al.* Antihypertensive effects and arterial hemodynamic alterations during angiotensin-converting enzyme inhibition. *J Hypertens* 1996;14:1139-46.
- 24 Asmar RG, London GM, Safar ME, for the REASON Project Coordinators and Investigators. Improvement in blood pressure, arterial stiffness and wave reflections with a very-low-dose Perindopril/Indapamide combination in hypertensive patient. *Hypertension* 2001;38:922-926.
- 25 Mahmud A, Feely J. Antihypertensive drugs and arterial stiffness. *Expert Rev Cardiovas Ther* 2003;1:65-78.
- 26 Morgan T, Lauri J, Bertram D, Anderson A. Effect of different antihypertensive drug classes on central aortic pressure. *Am J Hypertens* 2004;17:118-123.
- 27 Covic A, Goldsmith DJA, Gusbeth-Tatomir P, Covic M. Haemodialysis acutely improves endothelium-independent vasomotor function without significantly influencing the endothelium-mediated abnormal response to a B2-agonist. *Nephrol Dial Transplant* 2004;19:637-643.
- 28 Covic A, Goldsmith DJA, Panaghiu L, Covic M, Sedor J. Analysis of the effect of hemodialysis on peripheral and central arterial pressure waveforms. *Kidney Int* 2000;57:2634-2643.
- 29 Mardare N-G, Goldsmith DJA, Gusbeth-Tatomir P, Covic A. Intradialytic changes in reflective properties of the arterial system during a single hemodialysis session. *Hemodial Int* 2005;9:376-382
- 30 Vuurmans JLT, Boer WH, Bos WWW, Blankestijn PJ, Kooman HA. Contribution of volume overload and angiotensin II to the increased pulse wave velocity of haemodialysis patients. *J Am Soc Nephrol* 2002;13:177-183.
- 31 Covic A, Goldsmith DJA, Gusbeth-Tatomir P, Buhaescu I, Covic M. Successful renal transplantation decreases aortic stiffness and increases vascular reactivity in dialysis patients. *Transplantation* 2003;76:1573-1577.
- 32 Ferro CJ, Savage T, Pinder SJ, Tomson CRV. Central aortic pressure augmentation in stable renal transplant recipients. *Kidney Int* 2002;62:166-171.
- 33 Mustata S, Chan C, Lai V, Miller J. Impact of an exercise program on arterial stiffness and insulin resistance in hemodialysis patients. *J Am Soc Nephrol* 2004;5:2713–2718.
- 34 Covic A, Goldsmith DJA, Gusbeth-Tatomir P, Covic M. Haemodialysis acutely improves endothelium-independent vasomotor function without significantly influencing the endothelium-mediated abnormal response to a β 2-antagonist. *Nephrol Dial Transplant* 2004;19:637-643.